


MIRUS™ Ligating Clip is a sterile, medical grade titanium, non-absorbable, implantable device for general surgical procedures that require vessel or tissue ligation. **MIRUS™** Ligating Clip is intended to securely grip the vessel or tissue to provide complete occlusion, atraumatically.


MIRUS™ Ligating Clip are compatible with **MIRUS™** Ligating clip applicators, allowing the end user to ligate wide range of vessels and tissue structures.

Available as:

- Clip Size : Small, Medium, Medium-Large, Large
- Clips/Cartridges : 6
- Cartridges/Box : 20
- Clips/Box : 120
- Product Code : MLT100, MLT200, MLT300, MLT400



Meril Titanium Clips:



- Lateral & Transverse grooves allow secure fixation on the structure and increase resistance to dislodgement.
- Triangulated cross section maximises surface contact between clip & jaw.
- Precise tip to tip closure with firm vessel grip.

Ligation Clips

Titanium

Procedure Application

Blood vessel ligation for hemostasis

Appliers for Open Surgery:



Technical Specifications

Product	Product code	Specification	Length	Applicator size	Clips compatibility
Applicator-Open	FTOA-10015	OPEN APPLICATOR 100-15CM	15 cm	Small	MLT100
Applicator-Open	FTOA-10020	OPEN APPLICATOR 100-20CM	20 cm	Small	MLT100
Applicator-Open	FTOA-10028	OPEN APPLICATOR 100-28CM	28 cm	Small	MLT100
Applicator-Open	FTOA-20015	OPEN APPLICATOR 200-15CM	15 cm	Medium	MLT200
Applicator-Open	FTOA-20020	OPEN APPLICATOR 200-20CM	20 cm	Medium	MLT200
Applicator-Open	FTOA-20028	OPEN APPLICATOR 200-28CM	28 cm	Medium	MLT200
Applicator-Open	FTOA-30015	OPEN APPLICATOR 300-15CM	15 cm	Medium-Large	MLT300
Applicator-Open	FTOA-30020	OPEN APPLICATOR 300-20CM	20 cm	Medium-Large	MLT300
Applicator-Open	FTOA-30028	OPEN APPLICATOR 300-28CM	28 cm	Medium-Large	MLT300
Applicator-Open	FTOA-40015	OPEN APPLICATOR 400-15CM	15 cm	Large	MLT400
Applicator-Open	FTOA-40020	OPEN APPLICATOR 400-20CM	20 cm	Large	MLT400
Applicator-Open	FTOA-40028	OPEN APPLICATOR 400-28CM	28 cm	Large	MLT400

Appliers for Laparoscopic Surgery:



Technical Specifications

Product	Product code	Specification	Length	Applicator size	Compatible Clips	Shaft Diameter	Compatible Trocar
Applicator-Endo	FTEA-00300	LAPAROSCOPIC APPLICATOR 300	34 cm	Medium-Large	MLT300	10 mm	10 mm & above
Applicator-Endo	FTEA-00400	LAPAROSCOPIC APPLICATOR 400	34 cm	Large	MLT400	10 mm	10 mm & above