

Technical Specification

Stapler	Description	Open length	Close length	Cutting Diameter	Head Diameter	Rows of staples	Staple Quantity
MCS-21R3	Mirus Disp. Circular Stapler 3Row 21	4.5 mm	1.0-2.5 mm	12	21	3	21
MCS-24R3	Mirus Disp. Circular Stapler 3Row 24	4.5 mm	1.0-2.5 mm	15	24	3	24
MCS-25R3	Mirus Disp. Circular Stapler 3Row 25	4.5 mm	1.0-2.5 mm	16	25	3	33
MCS-26R3	Mirus Disp. Circular Stapler 3Row 26	4.5 mm	1.0-2.5 mm	17	26	3	33
MCS-29R3	Mirus Disp. Circular Stapler 3Row 29	4.5 mm	1.0-2.5 mm	20	29	3	36
MCS-32R3	Mirus Disp. Circular Stapler 3Row 32	4.5 mm	1.0-2.5 mm	22	32	3	42
MCS-24	Mirus Disp. Circular Stapler 24	4.5 mm	1.0-2.5 mm	15	24	2	18
MCS-25	Mirus Disp. Circular Stapler 25	4.5 mm	1.0-2.5 mm	16	25	2	20
MCS-26	Mirus Disp. Circular Stapler 26	4.5 mm	1.0-2.5 mm	17	26	2	20
MCS-29	Mirus Disp. Circular Stapler 29	4.5 mm	1.0-2.5 mm	20	29	2	24
MCS-32	Mirus Disp. Circular Stapler 32	4.5 mm	1.0-2.5 mm	22	32	2	30

Meril's Global Presence

Meril Tibbi Cihazlar

İmalat ve Ticaret A.Ş.
İçerenköy Mah. Çetinkaya Sok.
Prestij Plaza No:28, Kat:4 Ataşehir,
34752 İstanbul / Turkey
T : +90 216 641 44 24
F : +90 216 641 44 25

Meril Medical LLC.
Nauchnyi Proezd 19,
Moscow , Russia – 117 246.
Office - +7 495 772 7643

Meril UK Subsidiary
Crowe UK LLP, St. Brides
House, 10 Salisbury Square,
London EC4Y 8EH, UK

Meril GmbH.

Bornheimer Strasse 135-137,
D-53119 Bonn.
Germany.
T : +49 228 7100 4000
F : +49 228 7100 4001

Meril, inc.

2436 Emrick Boulevard,
Bethlehem, PA - 18020
T : (610)500-2080
F : (610)317-1672

Meril South America

Doc. Med LTDA Al. dos Tupiniquins,
1079 - Cep: 04077-003 - Moema.
Sao Paulo. Brazil.
T : +55 11 3624 5935
F : +55 11 3624 5936

EU Representative.
Obelis S.A.Bd, General Wahis 53,
1030, Brussels, Belgium.
T : +32 2 732 5954
F : +32 2 732 6003
E : mail@obelis.net

Meril SA Pty Ltd
102, 104, S101 and S102, Boulevard
West Office Park, 142 Western
Service Road, Erf 813
Woodmead Extension 17 Sandton,
Johannesburg – 2191 South Africa
T : +27 11 465-2049
F : +27 86 471 7941

Meril Endo Surgery Pvt. Ltd

Third Floor E1-E3, Meril Park,
Survey Number 135/2/B & 172/2
Muktanand Marg, Chala,
Vapi 396 191, Gujarat. India.
T : +91 260 240 8000
F : +91 260 305 2125
E : askinfo@merillife.com

Meril China Co. Ltd.

2301b,23f, Lixin Plaza,
no 90, South Hubin Road,
Xiamen, China.
T : 0086-592-5368505
F : 0086-592-5368519

Meril Bangladesh Pvt. Ltd.

Union Heights-01,
55-2, Bir Uttam Qazi
Nuruzzaman Sarak,
West Panthapath,
Dhaka-1205, Bangladesh
T : + 88 02 9122714/24

Meril Australia Pvt. Ltd.

Level 14, 60 Martin Place
Sydney NSW 2000, Australia
T : +61 28248 5800
F : +61 28248 5899

Meril Life Sciences Pvt. Ltd.

301, A - Wing, Business Square,
Chakala, Andheri Kurla Road,
Andheri East, Mumbai 400 093
T : + 91 22 39350700
F : + 91 22 39350777
E : askinfo@merillife.com


E askinfo@merillife.com
W www.merillife.com

Meril

Meril Endo-Surgery Pvt. Ltd.
E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala,
Vapi - 396 191, Gujarat, India

Meril

MIRUSTM
Circular Stapler




MIRUS™


Circular Stapler

MIRUS™ Circular Stapler is engineered for anastomosis in Gastrointestinal (GI) procedures. The design offers adequate lumen with a series of cutting diameters from 12mm to 22mm. With an adjustable staple height technology, it assures minimal tissue tension & controlled compression besides a reliable tactile & audible feedback while firing.

Meril is a global pioneer in 3 row stapling technology in endomechanical devices. MIRUS™ Circular Stapler comes from the same platform deploying an additional row of titanium staples thus creating better wound security & superior hemostasis. MIRUS™ Circular Stapler is also available in conventional 2 row model.



Ergonomic detachable anvil



Atraumatic dilating trocar tip



2 Row & 3 Row titanium staple lines ensuring superior hemostasis



White teflon cutting washer for audible feedback




Safety Lock prevents inadvertent firing of stapler



Wing Nut knob for controlled compression & tactile feedback



Optimal housing capacity to accommodate double doughnut



Tissue compression scale for adjustable staple height